



Optimizer®

Multiple Chamber Polyp Trap

The versatility of the Optimizer® Polyp Trap is highlighted by our independent, selectable specimen chambers, four direct bypass channels, and a convenient central vacuum lumen. A direct connection to wall suction makes for fast and easy fluid management.

Each trap is individually packaged for single-patient use to reduce risk of **cross-contamination**.



Features and Benefits

Versatile


- Four independent, selectable trap chambers
- Multiple direct bypass channels
- Visual markers for each chamber

Fast

- Convenient central vacuum lumen

Individually Packaged

- Single patient use reduces risk of cross-contamination



Direct connection to wall suction for easy fluid management



Sequenced trap chambers and direct bypass chambers

Ordering Information

Description	Qty	Catalog Number
Optimizer® Polyp Trap Multi Chamber Polyp Trap (Non-Sterile, Disposable)	20/Box	CAP132
Optimizer® Polyp Trap Multi Chamber Polyp Trap (Non-Sterile, Disposable)	120/Box	CAP132B

CONMED Corporation
11311 Concept Blvd, Largo, FL
33773

Toll Free: 1-866-4CONMED
International: 727-214-3000
www.CONMED.com

customerexperience@conmed.com
internationalorders@conmed.com