





Microbiology Brush

Prevents the risk of specimen contamination.



WANG
TBNA NEEDLE



PRECISOR
BRONCHO FORCEPS



PULMONARY
CYTOLOGY BRUSH



To learn more about this and other innovative products, call **800-448-6506**, visit **ConMed.com** or email **info.conmed.com**.

CONMED
MICROBIOLOGY
BRUSH



Sterility Protected Brush



Ring Handle


Sterile Specimens

ConMed has a full line of cytology brushes for your lab's sampling needs.

The ConMed 130 Microbiology Brush is designed to obtain microbiological specimens of the bronchial tree without contamination from the upper airway or bronchoscope.

Protected brushing allows the brush to pass through contaminated areas while maintaining sterility of the brush bristles.

A non-toxic, water soluble wax plug.



The Microbiology Brush features a non-toxic, water soluble wax plug, which prevents the risk of contamination of the brush prior to obtaining the specimen. After insertion, the plug is expelled and safely dissolved.

Microbiology Brush Ordering Information

Latex **FREE**

Catalog Number	Description	Brush Diameter	Length	Sheath Diameter	Unit
130	Microbiology Brush	1mm	110cm	1.8mm	5/box

Designed to obtain contamination free specimens, the brush is protected behind a wax plug which is expelled before sampling.



525 French Road
Utica, Ny 13502

Customer Service: 1800-448-6506
International Sales: +1(315) 797-8375
Fax: 1800-438-3051

Website: www.conmed.com
Email: info@mail.conmed.com

