

Bakterie-Virus/HME-filter



Bakterie/Virus/HME-filter



Bakterie/virusfilter

Produktinformation

- Komplett sortiment av bakterie-virus/HME-filter.
- Bakterieeffektivitet: 99.9999%, Viruseffektivitet 99.999%.

Material: SBC och PP



Orderinformation

| Artnr. | Beskrivning | Information | Avd. förp. |
|---------------|-------------------------------|--|------------|
| AL-08021.V001 | Bakterie-virus/HME filter vx. | Tidal volym: 150-1500ml Dead space: 55ml Bakterie: 99.9999% Virus: 99.999% Befuktning: 35,5 mgH ₂ O/L Vikt: 28,5gram | 50st |
| AL-08018.V001 | Bakterie-virus filter | Tidal volym: 150-1500ml Dead space: 33ml Bakterie: 99.9999% Virus: 99.999% Vikt: 19gram | 50st |
| AL-08020 | Bakterie-Virus HME-filter | Tidalvolym 150-1500ml Deadspace 53ml Bakterie: 99.9999% Virus: 99.999% Befuktning: 37,5 mgH ₂ O/l Vikt: 24 gram | 50st |

| | | | |
|---------------|--|---|------|
| AL-08031.V001 | Bakterie-Virus HME-filter pediatriskt | Tidal volym: 75.-300ml Dead space: 29ml Bakterie: 99.99% Virus: 99.99% Befuktning: Vikt: 21 gram | 75st |
|---------------|--|---|------|