



NOVO[®]

Everything from a single source

NOVO Product Catalogue 2023

Safe Solutions





Safe solutions for patient care

Innovation, creativity and technical knowledge are our strengths. This is how a product idea in 1987 became a company that is now one of the established and well-known manufacturers of medical fixations for intensive care.

Due to our proximity to the market and our claim to offer comprehensive and high-quality products for optimal clinical results, the expansion of our product range in the areas of ventilation, monitoring and patient safety was a logical step in 2018.

Today, we offer products for medical intensive care and perfectly fitting fixations from a single source. This is what we understand by genuine support for medical professionals in their daily clinical workflows.

Safe Solutions



Product development in-house

A deep technical understanding and a big heart for medical technology paired with a lot of creativity in the development of new product solutions, all this characterises our development department. In addition to the concerns of users and patients, details such as sustainable production and delivery strategies are also included in development processes.

Product and user trainings

Our consultants will show you what our products can do and how you can use them optimally for your applications. We know your requirements from practical experience and can therefore respond to your questions and wishes on site. Benefit from our expertise and make an appointment.

Our particular strengths

We are experts in what we do. We are constantly expanding our know-how and rely above all on the exchange with other experts. As a user, you know exactly what you need to optimise your daily work, how products should be designed, which individual adaptations would make sense for your patients. This is exactly where we are at your side as a partner. With professional advice, individual product optimisations and comprehensive training.

Individual product optimisation

Whether it's a small detail or a comprehensive product change: we can individually adapt selected products from our portfolio to your wishes and requirements. Contact us, in a personal discussion we will find the right solution together.



In-house production

Our medical fixations are still made by hand in-house. This means short distances and the greatest possible flexibility. From the processing of raw materials to packaging and shipping, quality assurance and fast response times are guaranteed.



Our Product Range

	Page
1. Respiratory Care	9
1.1 Airway Management	10
Endotracheal Tube Holders	10
Endotracheal Tubes & Accessories	18
Anaesthesia Masks	26
Laryngeal Masks	28
1.2 Breathing Circuits and Accessories	32
Dual-limb Breathing Circuits	36
Single-limb Breathing Circuits	38
Catheter Mounts	42
Rebreathing limbs with Breathing Bag	44
Breathing Circuits/Technical Data	46
1.3 Filtration & Humidification	48
Electrostatic Filters/HMEF	48
1.4 Tracheostomy	50
Tracheostomy Tube Holders	52
Tracheostoma Dressings	64
Wound Care	66
1.5 Inhalation & Aerosol Therapy	67
Nebulizers	67
1.6 Patient Positioning	72
Head Positioning Devices	72
2. Patient Safety	74
2.1 Decubitus Protection	75
Cushionings for Tubes	75
2.2 Phototherapy	76
Eye Protection	76
2.3 Patient Self-Protection	77
Hand and Foot Fixations	77
3. Infusion and Surgical Drain Management	80
3.1 Fixations	80
ECMO Fixations	82
Plaster Fixations	84
Joint Fixations	87
Fixations for Cables and Tubes	88
4. Monitoring	90
4.1 Blood Pressure Measurement	92
Fixations for Transducers	92
4.2 Cuff Pressure Measurement	93
Manometers	93
4.3 SPO ₂ Measurement	94
Fixations for Pulse Oximeter Sensors	94
5. Hygiene & Cleansing	96
5.1 Oral/Dental Care	96
Gels and Rinses	96
Oral Care Kits	96
Oral Swabs	96

Three colours - one variety





NOVO green

NOVO's own manufactured products - developed and manufactured in-house with the experience of over three decades.

NOVO green - the range of safe medical fixations.


A medical product is only as good as its fixation. With the help of the green colour dot you can quickly and easily find your desired product solution.



NOVO purple

NOVO branded quality products for respiratory care with a system


Everything from a single source - the product world with the NOVO purple colour dot offers the ideal complement for the care of your ventilated patients.



NOVO grey

Medical products to complement the NOVO range

The products marked in grey complete our green and purple range to form an overall portfolio.




Product Worlds



Product Groups



Product properties




1. Respiratory Care

In this section you will find solutions for all aspects of ventilation. The product range includes endotracheal tube holders, endotracheal tubes, filters, breathing circuits and comprehensive accessories.

1.1 Airway Management	10
1.2 Breathing Circuits & Accessories	32
1.3 Filtration & Humidification	48
1.4 Tracheostomy	50
1.5 Inhalation & Aerosol Therapy	67
1.6 Patient Positioning	72

1.1 Airway Management

Endotracheal Tube Holders



The NOVO range for endotracheal tube holders offers the ideal solution for all requirements. The soft and skin-friendly materials are gentle on the patient while positioning the tube securely. Find your preferred fixation technique to simplify your daily tasks.

One-piece, non-adhesive tube holders:

①



TuBo®-Clip



TuBo®-Clip S



TuBo®-Flex



Two-part, non-adhesive tube holders:

②



TuBo®-Fix



TuBo®-Grip




TuBo®-Guard Set



Multi-part adhesive tube holders:

③



TuBo®-Lock Set



TuBo®-Clip

One-piece endotracheal tube holder

Product properties

- One-piece tube holder with soft padding
- Fine foam structure on the skin side offers a good grip
- Appropriate for all types and sizes of endotracheal tubes
- Repositioning of the tube and refixation possible at any time
- Suitable for different head shapes, bearded patients and heavily sweating patients
- Suitable for MRI, non-adhesive
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft, laminated foam



REF	Article	Description	SU
20-01	TuBo®-Clip	Padding length 72 cm	100 pieces
20-05	TuBo®-Clip XL	Padding length 85 cm	100 pieces

TuBo®-Clip S

One-piece endotracheal tube holder

Product properties

- One-piece tube holder with soft padding
- Additional adhesive hook and loop strip for extra grip between tube and tube holder
- Fine foam structure on the skin side offers a good grip
- Appropriate for all types and sizes of endotracheal tubes
- Repositioning of the tube and refixation possible at any time
- Suitable for different head shapes, bearded patients and heavily sweating patients
- Suitable for MRI, non-adhesive
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft, laminated foam



REF	Article	Description	SU
20-02	TuBo®-Clip S	Padding length 72 cm	100 pieces
20-06	TuBo®-Clip S, XL	Padding length 85 cm	100 pieces

TuBo®-Flex

One-piece endotracheal tube holder

Product properties

- One-piece tube holder with soft padding
- With sewn-in taffeta strap for tube fixation
- Appropriate for all types and sizes of endotracheal tubes
- Easy and quick to use
- Can be removed without residue
- Suitable for different head shapes, bearded patients and heavily sweating patients
- Skin-friendly, soft materials
- Suitable for MRI, non-adhesive
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft, laminated foam and taffeta strap



REF	Article	Description	SU
20-30	TuBo®-Flex	Padding length 80 cm	100 pieces

TuBo®-Fix

Two-part endotracheal tube holder

Product properties

- Two-part tube holder consisting of middle part and neck padding
- Available in different versions and lengths
- Middle part, neck padding can be replaced individually if soiled
- Additional adhesive strips for extra security
- Repositioning of the tube and refixation possible at any time
- Can be cut individually
- Suitable for different head shapes, bearded patients and heavily sweating patients
- Skin-friendly, soft materials
- Suitable for MRI, non-adhesive
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft, laminated foam



REF	Article	Description	SU
20-10	TuBo®-Fix	Padding lengths 26 cm + 32.5 cm, 1 adhesive strip enclosed, hook and loop fastener	40 pieces
20-10K	TuBo®-Fix 1	Padding lengths 26 cm + 32.5 cm, 1 adhesive strip sewn in, hook and loop fastener	40 pieces
20-10-05	TuBo®-Fix 5	Padding lengths 28 cm + 32.5 cm, 1 adhesive strip sewn in, hook and loop fastener	40 pieces

Product accessories for TuBo®-Fix:

REF	Article	Description	SU
20-11	TuBo®-Fix middle part	Padding length 26 cm, 1 adhesive strip enclosed, component for TuBo®-Fix	80 pieces
20-11-05	TuBo®-Fix middle part 5	Padding length 28 cm, 1 adhesive strip sewn in, component for TuBo®-Fix 5	80 pieces
20-12	TuBo®-Fix middle part XL	Padding length 34 cm, 1 adhesive strip enclosed, extra long component for TuBo®-Fix	80 pieces
20-12-05	TuBo®-Fix middle part 5 XL	Padding length 34 cm, 1 adhesive strip sewn in, extra long component for TuBo®-Fix 5 XL	80 pieces
20-15	TuBo®-Fix neck padding	Padding length 32.5 cm, component for all TuBo®-Fix versions	80 pieces
20-16	TuBo®-Fix neck padding XL	Padding length 47 cm, extra long component for all TuBo®-Fix versions	80 pieces

TuBo®-Lock Set

Multi-part adhesive endotracheal tube holder

Product properties


- Multi-part tube holder consisting of middle part and 2 fixation plasters
- Available in different versions and lengths
- Components can be replaced individually
- Easy and quick to use
- Repositioning of the tube and refixation possible at any time
- Skin-friendly, soft materials
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft, laminated foam
- Plaster



REF	Article	Description	SU
25-10	TuBo®-Lock Set	Padding length 28 cm, adhesive strip affixed, fixation plasters	80 pieces
25-11	TuBo®-Lock Set XL	Padding length 28 cm, adhesive strip affixed, fixation plasters	80 pieces

Product accessories for TuBo®-Lock Set:

REF	Article	Description	SU
25-15	TuBo®-Lock middle part	Accessories for TuBo®-Lock Set, padding length 28 cm, adhesive strip affixed	80 pieces
25-16	TuBo®-Lock middle part XL	Accessories for TuBo®-Lock Set XL, padding length 36,5 cm, adhesive strip affixed	80 pieces
25-10-01	TuBo®-Lock fixation plasters	Fixation plasters for TuBo®-Lock Set	80 pieces

TuBo®-Grip

Two-part endotracheal tube holder

Product properties

- Two-part tube holder consisting of middle part with silicone-coated elastic strap on one side and neck padding
- Available in 2 different lengths
- Easy and quick to use
- Can be removed without leaving any residue
- Middle part, neck padding can be replaced individually if soiled
- Skin-friendly, soft materials
- Suitable for MRI, non-adhesive
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft, laminated foam and elastic strap



REF	Article	Description	SU
20-50	TuBo®-Grip	Middle part with silicone-coated elastic strap, padding lengths 33 cm + 32,5 cm, hook and loop fastener	60 pieces
20-51	TuBo®-Grip XL	Middle part with silicone-coated elastic strap, Padding lengths 33 cm + 47 cm, hook and loop fastener	60 pieces

Product accessories for TuBo®-Grip:

REF	Article	Description	SU
20-55	TuBo®-Grip middle part	Fixing accessories for TuBo®-Grip with elastic strap, padding length 33 cm, hook and loop fastener	60 pieces
20-52	TuBo®-Grip neck padding	Fixing accessories for TuBo®-Grip, padding length 32,5 cm	60 pieces
20-53	TuBo®-Grip neck padding XL	Fixing accessories for TuBo®-Grip, padding length 47 cm	60 pieces

TuBo®-Guard Set

Two-part endotracheal tube holder with bite guard

Product properties

- Two-part tube holder consisting of middle part with one-sided silicone coated elastic strap and neck padding as well as a compact bite guard to click on
- Available in different versions with colour-coded bite guard suitable for tube ID from 7.0 to 9.5 mm
- Easy and quick to use
- Can be removed without leaving any residue
- Middle part, neck padding can be replaced individually if soiled
- Skin-friendly, soft materials
- Suitable for MRI, non-adhesive
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft, laminated foam and bite guard



You can also find the Tubeguard bite guard individually on catalogue page 25.


REF	Article	Description	SU
20-50-7	TuBo®-Guard Set white	Padding lengths 33 cm + 32,5 cm, hook and loop fastener, bite guard white suitable for tube ID 7.0 mm	16 sets
20-50-75	TuBo®-Guard Set blue	Padding lengths 33 cm + 32,5 cm, hook and loop fastener, bite guard blue suitable for tube ID 7.5 mm	16 sets
20-50-8	TuBo®-Guard Set green	Padding lengths 33 cm + 32,5 cm, hook and loop fastener, bite guard green suitable for tube ID 8.0-8.5 mm	16 sets
20-50-9	TuBo®-Guard Set grey	Padding lengths 33 cm + 32,5 cm, hook and loop fastener, bite guard grey suitable for tube ID 9.0+9.5 mm	16 sets
20-51-7	TuBo®-Guard Set white XL	Padding lengths 33 cm + 47 cm, hook and loop fastener, bite guard white suitable for tube ID 7.0 mm	16 sets
20-51-75	TuBo®-Guard Set blue XL	Padding lengths 33 cm + 47 cm, hook and loop fastener, bite guard blue suitable for tube ID 7.5 mm	16 sets
20-51-8	TuBo®-Guard Set green XL	Padding lengths 33 cm + 47 cm, hook and loop fastener, bite guard green suitable for tube ID 8.0-8.5 mm	16 sets
20-51-9	TuBo®-Guard Set grey XL	Padding lengths 33 cm + 47 cm, hook and loop fastener, bite guard grey suitable for tube ID 9.0+9.5 mm	16 sets
20-51-X	TuBo®-Guard Set Mix XL	Padding lengths 33 cm + 47 cm, hook and loop fastener, bite guards white, blue, green and grey each 4 pieces	1 mixed carton

1.1 Airway Management

Endotracheal Tubes / Accessories

Anaesthesia Masks

AOVO®



Endotracheal Tube

Magill endotracheal tube with Murphy eye

Product properties

- Magill endotracheal tube
- Murphy eye, optionally with Hi-Lo cuff or without cuff
- Single packed
- Disposable product
- Sterile

Patient target group

- Neonates
- Small children
- Children
- Adolescents
- Adults

Material

- PVC
- Latex-free



REF with cuff	REF without cuff	Size	SU
IC-30 MC	IC-30 OC	3.0 mm ID	20 pieces
IC-35 MC	IC-35 OC	3.5 mm ID	20 pieces
IC-40 MC	IC-40 OC	4.0 mm ID	20 pieces
IC-45 MC	IC-45 OC	4.5 mm ID	20 pieces
IC-50 MC	IC-50 OC	5.0 mm ID	20 pieces
IC-55 MC	IC-55 OC	5.5 mm ID	20 pieces
IC-60 MC	IC-60 OC	6.0 mm ID	20 pieces
IC-65 MC	IC-65 OC	6.5 mm ID	20 pieces
IC-70 MC	-	7.0 mm ID	20 pieces
IC-75 MC	-	7.5 mm ID	20 pieces
IC-80 MC	-	8.0 mm ID	20 pieces
IC-85 MC	-	8.5 mm ID	20 pieces
IC-90 MC	-	9.0 mm ID	20 pieces

Endotracheal Tube

Spiral reinforced with stylet

Product properties

- Endotracheal tube with spiral reinforcement
- Murphy eye, optionally with Hi-Lo cuff or without cuff
- With stylet
- Single packed
- Disposable product
- Sterile

Patient target group

- Neonates
- Small children
- Children
- Adolescents
- Adults

Material

- PVC
- Metal
- Latex-free



REF with cuff	REF without cuff	Size	SU
ET1P30-R1MC	ET1P30-R1OC	3.0 mm ID	10 pieces
ET1P35-R1MC	ET1P35-R1OC	3.5 mm ID	10 pieces
ET1P40-R1MC	ET1P40-R1OC	4.0 mm ID	10 pieces
ET1P45-R1MC	ET1P45-R1OC	4.5 mm ID	10 pieces
ET1P50-R1MC	ET1P50-R1OC	5.0 mm ID	10 pieces
ET1P55-R1MC	ET1P55-R1OC	5.5 mm ID	10 pieces
ET1P60-R1MC	-	6.0 mm ID	10 pieces
ET1P65-R1MC	-	6.5 mm ID	10 pieces
ET1P70-R1MC	-	7.0 mm ID	10 pieces
ET1P75-R1MC	-	7.5 mm ID	10 pieces
ET1P80-R1MC	-	8.0 mm ID	10 pieces
ET1P85-R1MC	-	8.5 mm ID	10 pieces

Endotracheal Tube

Emergency Kit

Product properties

- Magill endotracheal tube set with stylet and 10 ml syringe
- Murphy-eye and Hi-Lo cuff
- Single packed / Disposable product
- Sterile


Patient target group

- Children
- Adolescents
- Adults

Material

- PVC
- Latex-free

REF	Size	SU
ET1P50-SL-SY	5.0 mm ID	10 pieces
ET1P55-SL-SY	5.5 mm ID	10 pieces
ET1P60-SL-SY	6.0 mm ID	10 pieces
ET1P65-SL-SY	6.5 mm ID	10 pieces
ET1P70-SL-SY	7.0 mm ID	10 pieces
ET1P75-SL-SY	7.5 mm ID	10 pieces
ET1P80-SL-SY	8.0 mm ID	10 pieces
ET1P85-SL-SY	8.5 mm ID	10 pieces



Endotracheal Tube

With suction line

Product properties

- Magill endotracheal tube with subglottic suction channel
- Murphy eye and Hi-Lo cuff
- Single packed / Disposable product
- Sterile


Patient target group

- Adults

Material

- PVC
- Latex-free

REF	Size	SU
217065 H	6.5 mm ID	10 pieces
217070H	7.0 mm ID	10 pieces
217075H	7.5 mm ID	10 pieces
217080H	8.0 mm ID	10 pieces
217085H	8.5 mm ID	10 pieces



Endotracheal Tube

Oral and nasal preformed

Product properties

- Oral or nasal preformed tube
- Murphy eye and Hi-Lo cuff
- Single packed
- Disposable product
- Sterile

Patient target group

- Neonates
- Small children
- Children
- Adolescents
- Adults

Material

- PVC
- Latex-free



Oral preformed:

REF with cuff	REF without cuff	Size	SU
212030H	212030U	3.0 mm ID	10 pieces
212035H	212035U	3.5 mm ID	10 pieces
212040H	212040U	4.0 mm ID	10 pieces
212045H	212045U	4.5 mm ID	10 pieces
212050H	212050U	5.0 mm ID	10 pieces
212055H	212055U	5.5 mm ID	10 pieces
212060H	212060U	6.0 mm ID	10 pieces
212065H	212065U	6.5 mm ID	10 pieces
212070H	-	7.0 mm ID	10 pieces
212075H	-	7.5 mm ID	10 pieces
212080H	-	8.0 mm ID	10 pieces
212085H	-	8.5 mm ID	10 pieces



Nasal preformed:

REF with cuff	REF without cuff	Size	SU
213030H	213030U	3.0 mm ID	10 pieces
213035H	213035U	3.5 mm ID	10 pieces
213040H	213040U	4.0 mm ID	10 pieces
213045H	213045U	4.5 mm ID	10 pieces
213050H	213050U	5.0 mm ID	10 pieces
213055H	213055U	5.5 mm ID	10 pieces
213060H	213060U	6.0 mm ID	10 pieces
213065H	-	6.5 mm ID	10 pieces
213070H	-	7.0 mm ID	10 pieces
213075H	-	7.5 mm ID	10 pieces
213080H	-	8.0 mm ID	10 pieces



Double Lumen Tube ○

For left or right side bronchial intubation

Product properties


- Double lumen tube left or right side
- With tracheal cuff and bronchial cuff (colour-coded)
- Accessories: suction catheter, connectors and Y-switch connector
- Single packed
- Disposable product
- Sterile

Patient target group

- Children
- Adolescents
- Adults

Material

- PVC
- Latex-free



REF left side	REF right side	Size	SU
231028	232028	28 Fr OD	1 piece
231032	232032	32 Fr OD	1 piece
231035	232035	35 Fr OD	1 piece
231037	232037	37 Fr OD	1 piece
231039	232039	39 Fr OD	1 piece
231041	232041	41 Fr OD	1 piece

Tubeguard ○

Bite guard

Product properties

- Colour-coded bite guard suitable for tube ID from 7.0 to 9.5 mm
- Handle for positioning the bite guard
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Thermoplastic elastomers and PP
- Latex-free



REF	Article	Description	SU
20-59-7	Tubeguard white	Bite guard white suitable for tube ID 7.0 mm	20 pieces
20-59-75	Tubeguard blue	Bite guard blue suitable for tube ID 7.5 mm	20 pieces
20-59-8	Tubeguard green	Bite guard green suitable for tube ID 8.0 - 8.5 mm	20 pieces
20-59-9	Tubeguard grey	Bite guard grey suitable for tube ID 9.0 - 9.5 mm	20 pieces

Shapable Intubation Stylet ○

For endotracheal tubes

Product properties


- Plastic shapable stylet in 3 different sizes
- Single packed
- Disposable product
- Sterile

Patient target group

- Children
- Adolescents
- Adults

Material

- Aluminium, PVC
- Latex-free



REF	Size	SU
SL1P06	6 Fr	40 pieces
SL1P10	10 Fr	40 pieces
SL1P14	14 Fr	40 pieces

Anaesthesia Mask ○

Product properties

- Anaesthesia mask with inflatable cuff and hook ring
- Available in different sizes
- Connection 22 mm ID
- Single packed / Disposable product
- Sterile

Patient target group


- Neonates
- Small children
- Children
- Adolescents
- Adults



Material

- PVC, PC
- Latex-free

REF with cuff	Size	Patient information	SU
IC-MA2P01-A	1	for neonates	100 pieces
IC-MA2P02-A	2	for small children	100 pieces
IC-MA2P03-A	3	for children and adults	100 pieces
IC-MA2P04A	4	for adults	100 pieces
IC-MA2P05-A	5	for adults	100 pieces
IC-MA2P06-A	6	for adults	100 pieces



1.1 Airway Management

Laryngeal Masks

The right choice for every requirement

In our Laryngeal Masks portfolio you will find the right product for every intubation. Our matrix shows you an overview of the different Laryngeal Masks variants.

Laryngeal Mask	Standard	Curved	Woodbrige	Inflatable Cuff	PVC
NOVO LM Standard PVC	✓			✓	✓
NOVO LM Curved PVC		✓		✓	✓
NOVO LM Woodbrige PVC			✓	✓	✓



NOVO LM Standard PVC

Disposable laryngeal mask PVC

Product properties

- PVC lumen
- Inflatable PVC cuff
- Connector 15 mm OD
- Single packed / Disposable product
- Sterile

Material

- PVC
- Latex-free



REF	Article	Details	Patient weight	SU
N10210010	NOVO LM Standard PVC 1	Size 1	≤ 5 kg	10 pieces
N10210015	NOVO LM Standard PVC 1.5	Size 1.5	5 – 10 kg	10 pieces
N10210020	NOVO LM Standard PVC 2	Size 2	10 – 20 kg	10 pieces
N10210025	NOVO LM Standard PVC 2.5	Size 2.5	20 – 30 kg	10 pieces
N10210030	NOVO LM Standard PVC 3	Size 3	30 – 50 kg	10 pieces
N10210040	NOVO LM Standard PVC 4	Size 4	50 – 70 kg	10 pieces
N10210050	NOVO LM Standard PVC 5	Size 5	> 70 kg	10 pieces

NOVO LM Curved PVC

Curved disposable laryngeal mask PVC

Product properties

- Curved PVC lumen
- Inflatable PVC cuff
- Connector 15 mm OD
- Single packed / Disposable product
- Sterile

Material

- PVC
- Latex-free



REF	Article	Details	Patient weight	SU
N10212010	NOVO LM Curved PVC 1	Size 1	≤ 5 kg	10 pieces
N10212015	NOVO LM Curved PVC 1.5	Size 1.5	5 – 10 kg	10 pieces
N10212020	NOVO LM Curved PVC 2	Size 2	10 – 20 kg	10 pieces
N10212025	NOVO LM Curved PVC 2.5	Size 2.5	20 – 30 kg	10 pieces
N10212030	NOVO LM Curved PVC 3	Size 3	30 – 50 kg	10 pieces
N10212040	NOVO LM Curved PVC 4	Size 4	50 – 70 kg	10 pieces
N10212050	NOVO LM Curved PVC 5	Size 5	70-100 kg	10 pieces
N10212060	NOVO LM Curved PVC 6	Size 6	> 100 kg	10 pieces

Comparison Standard - Curved



NOVO LM Woodbridge PVC

Disposable laryngeal mask PVC with spiral reinforcement

Product properties

- PVC lumen with spiral reinforcement
- Inflatable PVC cuff
- Connector 15 mm OD
- Single packed / Disposable product
- Sterile

Material

- PVC
- Metal
- Latex-free



REF	Article	Details	Patient weight	SU
N10214020	NOVO LM Woodbridge PVC 2	Size 2	10 – 20 kg	10 pieces
N10214025	NOVO LM Woodbridge PVC 2.5	Size 2.5	20 – 30 kg	10 pieces
N10214030	NOVO LM Woodbridge PVC 3	Size 3	30 – 50 kg	10 pieces
N10214040	NOVO LM Woodbridge PVC 4	Size 4	50 – 70 kg	10 pieces
N10214050	NOVO LM Woodbridge PVC 5	Size 5	> 70 kg	10 pieces

1.2 Breathing Circuits

Dual-limb Breathing Circuits

Single-limb Breathing Circuits

Catheter Mounts

Rebreathing limbs with Breathing Bag


The breathing circuit techniques at a glance

Extend extendable



- Lightweight
- Flexible in length
- Space saving storage

Smoothbore



- Very good drainage of the condensate
- Very good compliance due to smooth inner structure
- Kink resistant

Coaxial



- Single-limb system
- Inspiratory heating
- Very flexible
- Kink resistant

Bilumen 2-Chamber



- Single-limb system
- Lightweight
- Integrated gas sampling line (at „C“ equipment)
- Very good compliance
- Kink resistant

One equipment system - many options

With the NOVO breathing circuit matrix you will have everything at a glance. Find your right breathing circuit easily according to your personal requirements.

Equipment E
- Angled connector

Equipment CE
- Angled connector
- Gas sampling line

Equipment SE
- Angled connector
- Rebreathing limb with breathing bag

Equipment SE WT
- Angled connector
- Rebreathing limb with breathing bag
- Water trap

Equipment CSE
- Angled connector
- Gas sampling line
- Rebreathing limb with breathing bag

Equipment CSE WT
- Angled connector
- Gas sampling line
- Rebreathing limb with breathing bag
- Water trap

Equipment HE
- Angled connector
- Angled connector
- HMEF
- Catheter mount

Angled connector Gas sampling line Filter/HMEF Water trap

Catheter mount Rebreathing limb with breathing bag

Fields of application		Dual-Limb Breathing Circuits				Single-Limb Breathing Circuits					
		Extend Extendable Limb		Smoothbore Smoothbore Limb		Coaxial Tube-in-Tube System			Bilumen Breathing Circuit Two Chamber		
Equipment	Fields of application	180 cm	240 cm	180 cm	240 cm	180 cm	240 cm	300 cm	180 cm	240 cm	300 cm
E	Anesthesia	N20418010	N20424010	N20318030	N20324010	N20118010	N20124010		IC-QT 180 E	IC-QT 240 E	IC-QT 300 E
CE	Intensive Care								IC-QT 180 CE	IC-QT 240 CE	IC-QT 300 CE
SE			N20424030			N20118030	N20124030		IC-QT 180 SE		
SE WT					N20324035						
CSE						N20118040	N20124040	N20130040	IC-QT 180 CSE	IC-QT 240 CSE	IC-QT 300 CSE
CSE WT									IC-QT 180 CSE WT	IC-QT 240 CSE WT	
HE						N20118017					
E Kids		N20418011	N20424011								
SE Kids			N20424031								

● expiratory ● inspiratory

NOVO BC Extend

Extendable dual-limb breathing circuit

E/SE

NEW

Product properties

- Angled connector
- Norm connectors 22 mm OD/15 mm ID
- Single packed / Disposable product

Equipment E

- Extendable up to 1,8 m / 2,4 m

Equipment SE

- Extendable up to 2,4 m
- Smoothbore rebreathing limb 0,9 m
- Breathing Bag 2,0 l



REF	Article	Details	SU
N20418010	NOVO BC Extend 180E	Length 180 cm	20 pieces
N20424010	NOVO BC Extend 240E	Length 240 cm	20 pieces
N20424030	NOVO BC Extend 240SE	Length 240 cm	20 pieces

NOVO BC Extend Kids

Extendable dual-limb breathing circuit for children

E/SE

NEW

Product properties

- Angled connector
- Norm connectors 22 mm OD/15 mm ID
- Single packed / Disposable product

Equipment E

- Extendable up to 1,8 m / 2,4 m

Equipment SE

- Extendable up to 2,4 m
- Smoothbore rebreathing limb 0,9 m
- Breathing Bag 1,0 l



REF	Article	Details	SU
N20418011	NOVO BC Extend Kids 180E	Length 180 cm	20 pieces
N20424011	NOVO BC Extend Kids 240E	Length 240 cm	20 pieces
N20424031	NOVO BC Extend Kids 240SE	Length 240 cm	20 pieces

NOVO BC Smoothbore

Smoothbore breathing circuit

E/SE WT

NEW

Product properties

- Angled connector
- Norm connectors 22 mm OD/15 mm ID
- Single packed / Disposable product

Equipment E

- Length 1,8 m / 2,4 m

Equipment SE WT

- Length 2,4 m
- Smoothbore rebreathing limb 0,9 m
- Breathing Bag 2,0 l
- Double water trap



Equipment SE WT:
Double water trap and
rebreathing limb and breathing bag

REF	Article	Details	SU
N20318030	NOVO BC Smoothbore 180E	Length 180 cm	40 pieces
N20324010	NOVO BC Smoothbore 240E	Length 240 cm	35 pieces
N20324035	NOVO BC Smoothbore 240SE WT	Length 240 cm	25 pieces

NOVO BC QTube™ Bilumen

E/SE/CE

Two-chamber breathing circuit

Product properties

- Angled connector
- Norm connectors 22 mm OD/15 mm ID
- Single packed / Disposable product

Equipment E



- Length 1,8 m / 2,4 m / 3,0 m

Equipment SE



- Length 1,8 m
- Corrugated rebreathing limb 1,2 m
- Breathing Bag 2,0 l

Equipment CE



- Length 1,8 m / 2,4 m / 3,0 m
- Integrated gas sampling line
- External gas sampling line 1,2 m"



Integrated gas sampling line



Additional external gas sampling line

REF	Article	Details	SU
IC-QT180 E	NOVO BC QTube™ Bilumen 180E	Length 180 cm	30 pieces
IC-QT 240 E	NOVO BC QTube™ Bilumen 240E	Length 240 cm	30 pieces
IC-QT 300 E	NOVO BC QTube™ Bilumen 300E	Length 300 cm	30 pieces
IC-QT 180 SE	NOVO BC QTube™ Bilumen 180SE	Length 180 cm	30 pieces
IC-QT 180 CE	NOVO BC QTube™ Bilumen 180CE	Length 180 cm	30 pieces
IC-QT 240 CE	NOVO BC QTube™ Bilumen 240CE	Length 240 cm	30 pieces
IC-QT 300 CE	NOVO BC QTube™ Bilumen 300CE	Length 300 cm	30 pieces

NOVO BC QTube™ Bilumen

CSE/CSE WT

Two-chamber breathing circuit

Product properties

- Angled connector
- Norm connectors 22 mm OD/15 mm ID
- Single packed / Disposable product

Equipment CSE



- Length 1,8 m / 2,4 m / 3,0 m
- Integrated gas sampling line
- External gas sampling line 1,2 m
- Corrugated rebreathing limb 1,2 m
- Breathing Bag 2,0 l

Equipment CSE WT



- Length 1,8 m / 2,4 m
- Integrated gas sampling line
- External gas sampling line 1,2 m
- Corrugated rebreathing limb 1,2 m
- Breathing Bag 2,0 l
- Double water trap



Double water trap (Equipment WT)



Integrated gas sampling line



Additional external gas sampling line

REF	Article	Details	SU
IC-QT180 CSE	NOVO BC QTube™ Bilumen 180CSE	Length 180 cm	30 pieces
IC-QT 240 CSE	NOVO BC QTube™ Bilumen 240CSE	Length 240 cm	25 pieces
IC-QT 300 CSE	NOVO BC QTube™ Bilumen 300CSE	Length 300 cm	20 pieces
IC-QT 180 CSE WT	NOVO BC QTube™ Bilumen 180CSE WT	Length 180 cm	20 pieces
IC-QT 240 CSE WT	NOVO BC QTube™ Bilumen 240CSE WT	Length 240 cm	20 pieces

NOVO BC Coax

Coaxial breathing circuit

E/SE/CSE

NEW

Product properties

- Angled connector
- Norm connectors 22 mm OD/15 mm ID
- Single packed / Disposable product

Equipment E



- Length 1,8 m / 2,4 m

Equipment SE



- Length 1,8 m / 2,4 m
- Rebreathing limb 0,9 m
- Breathing bag 2,0 l

Equipment CSE



- Length 1,8 m / 2,4 m / 3,0 m
- External gas sampling line 1,2 m
- Rebreathing limb 0,9 m
- Breathing bag 2,0 l



REF	Article	Details	SU
N20118010	NOVO BC Coax 180E	Length 180 cm	20 pieces
N20124010	NOVO BC Coax 240E	Length 240 cm	30 pieces
N20118030	NOVO BC Coax 180SE	Length 180 cm	30 pieces
N20124030	NOVO BC Coax 240SE	Length 240 cm	25 pieces
N20118040	NOVO BC Coax 180CSE	Length 180 cm	30 pieces
N20124040	NOVO BC Coax 240CSE	Length 240 cm	25 pieces
N20130040	NOVO BC Coax 300CSE	Length 300 cm	25 pieces

NOVO BC Coax Set HE

Coaxial breathing circuit

HE E

NEW


Product properties

- Angled connector
- Norm connectors 22 mm OD/15 mm ID
- Single packed / Disposable product

Equipment HE



- Length 1,8 m
- Smoothbore catheter mount 15 cm
- HMEF



REF	Article	Details	SU
N20118017	NOVO BC Coax Set HE 180E	Length 180 cm	30 pieces



NOVO CM Extend S

Extendable Catheter mount

NEW

Extendable catheter mount for use in children, adolescents and adults.

Product properties

Version with angled connector

- With double swivel connector, bronchoscopy and suction port
- Adapter 22 mm OD/15 mm ID - 22 mm ID
- Length 15 cm
- Single packed / Disposable product
- Sterile



Version with straight connector

- Adapter 22 mm OD/15 mm ID - 22 mm ID
- Length 15 cm
- Single packed / Disposable product
- Sterile

REF	Article	Details	SU
N22220410	NOVO CM 15 Extend S	Length 15 cm / Double swivel connector	50 pieces
N22220411	NOVO CM 15 Extend Straight S	Length 15 cm straight	50 pieces

NOVO CM Connector Double Swivel S

Double swivel connector

NEW

Product properties

- Double swivel connector, bronchoscopy and suction port
- Adapter 22 mm OD/15 mm ID - 15 mm OD
- Single packed / Disposable product
- Sterile



REF	Article	SU
N22220412	NOVO CM Connector Double Swivel S	75 pieces

NOVO CM Smoothbore S


Smoothbore catheter mount

NEW

Smoothbore catheter mount for use in children, adolescents and adults.

Product properties

- With double swivel connector, bronchoscopy and suction port
- Adapter 22 mm OD/15 mm ID - 22 mm ID
- Optionally available in 15 cm or 20 cm length
- Single packed / Disposable product
- Sterile



REF	Article	Details	SU
N22220310	NOVO CM 15 Smoothbore S	Length 15 cm / Double swivel connector	50 pieces
N2222020	NOVO CM 20 Smoothbore S	Length 20 cm / Double swivel connector	50 pieces

NOVO BA Breathing Bag

Smoothbore rebreathing limb with breathing bag

NEW

Product properties

- Breathing bag 1,0 l or 2,0 l
- Smoothbore limb, length 1 m
- Single packed / Disposable product



REF	Article	Details	SU
N22010001	NOVO BA Breathing Bag 1/100	Volume 1,0 l / Limb length 100 cm	50 pieces
N22010002	NOVO BA Breathing Bag 2/100	Volume 2,0 l / Limb length 100 cm	50 pieces


NOVO BA Breathing Bag Corrugated

Corrugated rebreathing limb with breathing bag


NEW

Product properties

- Breathing bag 0,5 l / 1,0 l oder 2,0 l
- Corrugated limb, length 1 m
- Single packed / Disposable product



REF	Article	Details	SU
N22010105	NOVO BA Breathing Bag Corrugated 0.5/100	Volume 0,5 l / Limb length 100 cm	50 pieces
N22010101	NOVO BA Breathing Bag Corrugated 1/100	Volume 1,0 l / Limb length 100 cm	50 pieces
N22010102	NOVO BA Breathing Bag Corrugated 2/100	Volume 2,0 l / Limb length 100 cm	50 pieces



Breathing Circuits/Technical Data

NOVO BC Extend Extendable dual-limb breathing circuit

Equipment E

NOVO BC Extend 180E REF N20418010
NOVO BC Extend 240E REF N20424010



Equipment SE

NOVO BC Extend 240SE REF N20424030



NOVO BC Extend Kids Extendable dual-limb breathing circuit


Equipment E

NOVO BC Extend Kids 180E REF N20418011
NOVO BC Extend Kids 240E REF N20424011



Equipment SE

NOVO BC Extend Kids 240SE REF N20424031



NOVO BC Smoothbore Smoothbore breathing circuit


Equipment E

NOVO BC Smoothbore 180E REF N20318030
NOVO BC Smoothbore 240E REF N20324010



Equipment SE WT


NOVO BC Smoothbore 240SE WT REF N20324035



NOVO BC Coax / NOVO BC Coax Set Coaxial breathing circuit


Equipment E

NOVO BC Coax 180E REF N20118010
NOVO BC Coax 240E REF N20124010




Equipment SE

NOVO BC Coax 180SE REF N20118030
NOVO BC Coax 240SE REF N20124030




Equipment CSE

NOVO BC Coax 180CSE REF N20118040
NOVO BC Coax 240CSE REF N20124040
NOVO BC Coax 300CSE REF N20130040



Equipment HE

NOVO BC Coax Se 180HE REF N20118017




Breathing Circuits/Technical Data

NOVO BC QTube™ Bilumen Two-chamber breathing circuit


Equipment E

NOVO BC QTube™ Bilumen 180E REF IC-QT 180 E
NOVO BC QTube™ Bilumen 240E REF IC-QT 240 E
NOVO BC QTube™ Bilumen 300E REF IC-QT 300 E




Equipment CE

NOVO BC QTube™ Bilumen 180CE REF IC-QT 180 CE
NOVO BC QTube™ Bilumen 240CE REF IC-QT 240 CE
NOVO BC QTube™ Bilumen 300CE REF IC-QT 300 CE




Equipment SE

NOVO BC QTube™ Bilumen 180SE REF IC-QT 180 SE




Equipment CSE

NOVO BC QTube™ Bilumen 180CSE REF IC-QT 180 CSE
NOVO BC QTube™ Bilumen 240CSE REF IC-QT 240 CSE
NOVO BC QTube™ Bilumen 300CSE REF IC-QT 300 CSE



Equipment CSE WT

NOVO BC QTube™ Bilumen 180CSE WT REF IC-QT 180 CSE WT
NOVO BC QTube™ Bilumen 240CSE WT REF IC-QT 240 CSE WT



NOVO BA Breathing Bag Smoothbore rebreathing limb with breathing bag

NOVO BA Breathing Bag 1/100 REF N22010001
NOVO BA Breathing Bag 2/100 REF N22010002




NOVO BA Breathing Bag Corrugated Corrugated rebreathing limb with breathing bag


NOVO BA Breathing Bag Corrugated 0,5/100 REF N22010105
NOVO BA Breathing Bag Corrugated 1/100 REF N22010101
NOVO BA Breathing Bag Corrugated 2/100 REF N22010102



NOVO CM Extend / Straight S Catheter mount




NOVO CM 15 Extend S REF N22220410



NOVO CM 15 Extend Straight S REF N22220411

NOVO CM Smoothbore S Catheter mount

NOVO CM 15 Smoothbore REF N22220310
NOVO CM 20 Smoothbore REF N2222020





NOVO FI Breathing Filter S

Bacterial and viral filter

Product properties

- With Luer Port
- Single packed / Disposable product
- Sterile

Tidal Volume:	150 – 1500 ml
Resistance to Flow at:	30 l/min: ≤ 0.8 cmH ₂ O 60 l/min: ≤ 1.7 cmH ₂ O 90 l/min: ≤ 2.7 cmH ₂ O
Dead Space:	33 ml
Mode of Filtration:	Elektrostatic
Bacterial Filtration Efficiency:	99,999 %
Viral Filtration Efficiency:	99,9 %
Connectors:	22mmM 15 / 22 mm MF

REF	Article	Details	SU
N16410010	NOVO FI Breathing Filter S	With Luer Port	50 pieces



NEW

NOVO FI HMEF S

HME filter

Product properties

- With Luer Port
- Single packed / Disposable product
- Sterile

Tidal Volume:	150 – 1500 ml
Resistance to Flow at:	30 l/min: ≤ 1.0 cmH ₂ O 60 l/min: ≤ 2.2 cmH ₂ O 90 l/min: ≤ 3.6 cmH ₂ O
Dead Space:	53 ml
Mode of Filtration:	Elektrostatisch
Bacterial Filtration Efficiency:	99,999 %
Viral Filtration Efficiency:	99,9 %
Humidification Efficiency:	37,3 mg/l H ₂ O (500 ml tid. vol).
Connectors:	22mmM 15 / 22 mm MF

REF	Article	Details	SU
N16110012	NOVO FI HMEF S	With Luer Port	45 pieces



NEW

1.3 Filtration & Humidification

Electrostatic Filters/HMEF

1.4 Tracheostomy

Tracheostomy Tube Holders

Tracheostoma Dressings

Wound Care



Individual Product Optimisations

We can individually customise selected products from our portfolio to your wishes and requirements. You can recognise them by our „Optimisation“ icon.

Possible optimisations:

- Length adjustments
- Width of the elastic straps
- Partition of the straps


Please feel free to contact us, together we will find the optimal solution for you.

How to determine the correct length of the tracheostomy tube holder


Measure the circumference of the patient's neck up to the tube holder slots of the cannula shield.

Compare this measurement with the indicated length ranges of the tube holders.

If your measured value is within the specified range, you have the optimum tube holder size for your patient.



Reference points for the length specifications of the NOVO tracheostomy tube holders



1.4 Tracheostomy

Tracheostomy Tube Holders



One-piece tube holders:

1




NOVO Edera®
 ○ 1 ■■■■
 ♻️ ⚙️ ⚙️ ⚙️



BePa-Clip®
 ○ 1 ♻️
 ⚙️ ⚙️



Trachea-Fix®
 ○ 1 ■■■■
 ♻️ ⚙️



Trach-Basic
 ○ 1 ■■■■ ⚙️



Sim-Plex
 ○ 1 ⚙️

One-piece tube holders for children:


1 🐻



NOVO Edera® Kids
 ○ 1 ■■■■ ♻️
 ⚙️ ⚙️ 🐻



BePa-KidsClip®
 ○ 1 ♻️
 ⚙️ ⚙️ 🐻



Trachea-Fix® Kids
 ○ 1 ■■■■ ♻️
 ♻️ ⚙️ 🐻

Two-part tube holders:


2



NOVO Edera® Duo
 ○ 2 ■■■■ ♻️
 ♻️ ⚙️ ⚙️ ⚙️



Duo-Clip®
 ○ 2 ■■■■ ♻️
 ♻️ ⚙️ ⚙️



Duo-Fix
 ○ 2 ■■■■ ♻️
 ♻️ ⚙️

Two-part tube holders for children:

2 🐻



NOVO Edera® Duo Kids
 ○ 2 ■■■■ ♻️ ⚙️ ⚙️ 🐻



Duo-KidsClip
 ○ 2 ■■■■ ♻️ ⚙️ 🐻

NOVO Edera®



One-piece tracheostomy tube holder

Product properties

- Tracheostomy tube holder with velour lamination, 2,5 cm padding width
- Body side with skin-friendly cotton lamination
- Available with hook or hook and loop fastener
- Centric positioning of the padding due to length adjustment on both sides
- Easy handling
- Soft materials
- Suitable for MRI
- Also available with individual optimisation
- Disposable product

Patient target group

- Children
- Adolescents
- Adults

Material

- Soft foam with velour and cotton lamination
- Elastic strap



REF 1 piece	REF 30 pieces	Article	Description
N15010018	N15030018	NOVO Edera® 18/H	Length range 20,5 – 39,0 cm, 18 cm padding, hook fastener
N15010025	N15030025	NOVO Edera® 25/H	Length range 28,0 – 47,0 cm, 26 cm padding, hook fastener
N15010036	N15030036	NOVO Edera® 36/H	Length range 38,0 – 57,0 cm 36 cm padding, hook fastener
N15010118	N15030118	NOVO Edera® 18/K	Length range 22,5 – 41,0 cm, 18 cm padding, hook and loop fastener
N15010125	N15030125	NOVO Edera® 25/K	Length range 29,5 – 49,0 cm, 26 cm padding, hook and loop fastener
N15010136	N15030136	NOVO Edera® 36/K	Length range 40,0 – 59,0 cm, 36 cm padding, hook and loop fastener

NOVO Edera® tube holders with optimisations according to your special requirements

REF 1 piece	REF 30 pieces	Article	Description
N15110001	N15110030	NOVO Edera® Opti H	Optimisation tube holder with hook fastener
N15120001	N15120030	NOVO Edera® Opti K	Optimisation tube holder with hook and loop fastener

NOVO Edera® Kids



One-piece tracheostomy tube holder for children

Product properties

- Tracheostomy tube holder with velour lamination, 2,5 cm padding width
- Body side with skin-friendly cotton lamination
- Available with hook and loop fastener
- Centric positioning of the padding due to length adjustment on both sides
- Easy handling
- Soft materials
- Suitable for MRI
- Also available with individual optimisation
- Disposable product

Patient target group

- Premature babies
- Neonates
- Infants
- Small children

Material

- Soft foam with velour and cotton lamination
- Elastic strap



REF 1 piece	REF 30 pieces	Article	Description
N15010107	N15030107	NOVO Edera® Kids 7/K	Length range 12,0 – 21,0 cm, 7,5 cm padding, hook and loop fastener
N15010111	N15030111	NOVO Edera® Kids 11/K	Length range 15,0 – 24,0 cm, 11 cm padding, hook and loop fastener
N15010116	N15030116	NOVO Edera® Kids 16,5/K	Length range 20,0 – 30,0 cm, 16,5 cm padding, hook and loop fastener

NOVO Edera® tube holders with optimisations according to your special requirements

REF 1 piece	REF 30 pieces	Article	Description
N15110001	N15110030	NOVO Edera® Opti H	Optimisation tube holder with hook fastener
N15120001	N15120030	NOVO Edera® Opti K	Optimisation tube holder with hook and loop fastener

NOVO Edera® Duo

Two-part tracheostomy tube holder

Product properties


- Tracheostomy tube holder with velour lamination, 2,5 cm padding width
- Body side with skin-friendly cotton lamination
- Available with hook or hook and loop fastener
- Centric positioning of the padding due to length adjustment on both sides
- Easy handling
- Soft materials
- Suitable for MRI
- Also available with individual optimisation
- Disposable product

Patient target group

- Children
- Adolescents
- Adults

Material

- Soft foam with velour and cotton lamination
- Elastic strap



NOVO Edera® Duo in detail:

Two-part tube holder for an easy and comfortable application.

REF 1 piece	REF 30 pieces	Article	Description
N15020018	N15040018	NOVO Edera® Duo 18/H	Length range 20,5 – 39,0 cm, 18 cm padding, hook fastener
N15020025	N15040025	NOVO Edera® Duo 25/H	Length range 28,0 – 47,0 cm, 26 cm padding, hook fastener
N15020036	N15040036	NOVO Edera® Duo 36/H	Length range 38,0 – 57,0 cm, 36 cm padding, hook fastener
N15020118	N15040118	NOVO Edera® Duo 18/K	Length range 22,5 – 41,0 cm, 18 cm padding, hook and loop fastener
N15020125	N15040125	NOVO Edera® Duo 25/K	Length range 29,5 – 49,0 cm, 26 cm padding, hook and loop fastener
N15020136	N15040136	NOVO Edera® Duo 36/K	Length range 40,0 – 59,0 cm, 36 cm padding, hook and loop fastener

NOVO Edera® tube holders with optimisations according to your special requirements

REF 1 piece	REF 30 pieces	Article	Description
N15110001	N15110030	NOVO Edera® Opti H	Optimisation tube holder with hook fastener
N15120001	N15120030	NOVO Edera® Opti K	Optimisation tube holder with hook and loop fastener

NOVO Edera® Duo Kids

Two-part tracheostomy tube holder for children

Product properties


- Tracheostomy tube holder with velour lamination, 2,5 cm padding width
- Body side with skin-friendly cotton lamination
- Available with hook and loop fastener
- Centric positioning of the padding due to length adjustment on both sides
- Easy handling
- Soft materials
- Suitable for MRI
- Also available with individual optimisation
- Disposable product

Patient target group

- Infants
- Small children

Material

- Soft foam with velour and cotton lamination
- Elastic strap



NOVO Edera® Duo Kids in detail:

Two-part tube holder for an easy and comfortable application.

REF 1 piece	REF 30 pieces	Article	Description
N15020116	N15040116	NOVO Edera® Duo Kids 16,5/K	Length range 20,0 – 30,0 cm, 16,5 cm padding, hook and loop fastener

NOVO Edera® tube holders with optimisations according to your special requirements

REF 1 piece	REF 30 pieces	Article	Description
N15110001	N15110030	NOVO Edera® Opti H	Optimisation tube holder with hook fastener
N15120001	N15120030	NOVO Edera® Opti K	Optimisation tube holder with hook and loop fastener

BePa-Clip®

One-piece tracheostomy tube holder



Product properties

- Tracheostomy tube holder with 3 cm wide and breathable padding
- Available with hook or hook and loop fastener
- Centric positioning of the padding due to length adjustment on both sides
- Easy handling
- Skin-friendly, soft materials
- Suitable for MRI
- Disposable product

Patient target group

- Children
- Adolescents
- Adults

Material

- Soft foam with skin-friendly coating
- Elastic strap
- Latex-free



REF 30 pieces*	REF 60 pieces	Article	Description
100 444	200 444	BePa-Clip® STANDARD 1	Length range 29,0 – 34,8 cm, 26 cm padding, hook fastener
-	200 443	BePa-Clip® STANDARD 2	Length range 36,0 – 42,0 cm, 26 cm padding, hook fastener
100 301	200 301	BePa-Clip® VARIO 18/H	Length range 21,5 – 38,0 cm, 18 cm padding, hook fastener
100 302	200 302	BePa-Clip® VARIO 18/K	Length range 23,0 – 38,0 cm, 18 cm padding, hook and loop fastener
100 501	200 501	BePa-Clip® VARIO 25/H	Length range 29,0 – 47,0 cm, 25 cm padding, hook fastener
100 502	200 502	BePa-Clip® VARIO 25/K	Length range 29,5 – 48,8 cm 25 cm padding, hook and loop fastener
100 601	200 601	BePa-Clip® VARIO 30/H	Length range 33,5 – 51,0 cm, 30 cm padding, hook fastener
100 602	200 602	BePa-Clip® VARIO 30/K	Length range 35,0 – 51,0 cm, 30 cm padding, hook and loop fastener
100 701	200 701	BePa-Clip® VARIO 36/H	Length range 39,0 – 56,0 cm, 36 cm padding, hook fastener
100 702	200 702	BePa-Clip® VARIO 36/K	Length range 41,0 – 57,0 cm, 36 cm padding, hook and loop fastener

* single packed

BePa-KidsClip®

One-piece tracheostomy tube holder for children



Product properties

- Tracheostomy tube holder with 3 cm wide and breathable padding
- Available with hook or hook and loop fastener
- Centric positioning of the padding due to length adjustment on both sides
- Easy handling
- Skin-friendly, soft materials
- Suitable for MRI
- Disposable product

Patient target group

- Premature babies
- Neonates
- Infants
- Small children

Material

- Soft foam with skin-friendly coating
- Elastic strap
- Latex-free



REF 30 pieces*	REF 60 pieces	Article	Description
-	SA 7/K	BePa-KidsClip® 7/K	Length range 12,0 – 16,0 cm, 7 cm padding, hook and loop fastener
-	SA 9/K	BePa-KidsClip® 9/K	Length range 13,5 – 20,5 cm, 9 cm padding, hook and loop fastener
100 201	200 201	Bepa-KidsClip® 11/H	Length range 14,7 – 24,5 cm, 11 cm padding, hook fastener
100 202	200 202	BePa-KidsClip® 11/K	Length range 16,5 – 24,0 cm, 11 cm padding, hook and loop fastener
100 211	200 211	BePa-KidsClip® 13/H	Length range 15,7 – 24,5 cm, 13 cm padding, hook fastener
100 212	200 212	BePa-KidsClip® 13/K	Length range 17,0 – 24,2 cm, 13 cm padding, hook and loop fastener
100 221	200 221	BePa-KidsClip® 15/H	Length range 18,0 – 27,5 cm, 15 cm padding, hook fastener
100 222	200 222	BePa-KidsClip® 15/K	Length range 19,5 – 27,0 cm, 15 cm padding, hook and loop fastener

* single packed

Trachea-Fix®

One-piece tracheostomy tube holder

Product properties

- Tracheostomy tube holder with soft 3 cm wide padding
- Strap length adjustable from 21-42 cm
- Easy and quick to use
- Suitable for different neck circumferences
- Hook and loop fastener with rounded hook and loop strap for quick and easy closing
- Skin-friendly, soft and breathable materials
- Dimensionally stable - even when worn for longer periods of time
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft foam with skin-friendly lamination
- Latex-free



REF	Article	Description	SU
200 442	Trachea-Fix®	Length range 21 - 42 cm, strap length 45 cm, strap width 3 cm	60 pieces
100 442	Trachea-Fix®	Length range 21 - 42 cm, strap length 45 cm, strap width 3 cm, single packed	30 pieces

Trachea-Fix® Kids

One-piece tracheostomy tube holder for children

Product properties

- Tracheostomy tube holder with soft 2 cm wide padding
- Variably adjustable
- Easy and quick to use
- Suitable for different neck circumferences
- Hook and loop fastener with rounded hook and loop strap for quick and easy closing
- Skin-friendly, soft and breathable materials
- Dimensionally stable - even when worn for longer periods of time
- Suitable for MRI
- Disposable product

Patient target group

- Small children
- Children

Material

- Soft foam with skin-friendly lamination
- Latex-free



REF	Article	Description	SU
200 542	Trachea-Fix® Kids	Strap length 25,5 cm, strap width 2 cm	60 pieces

Trach-Basic

One-piece tracheostomy tube holder

Product properties

- Tracheostomy tube holder with soft 3,5 cm wide padding
- Strap length adjustable from 17-46 cm
- Hook and loop fastener
- Easy and quick to use
- Suitable for different neck circumferences
- Skin-friendly, soft and breathable materials
- Inconspicuously to wear
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft foam with skin-friendly lamination
- Latex-free



REF	Article	Description	SU
100 441	Trach-Basic	Length range 17 - 46 cm, Strap length 46,5 cm, strap width 3,5 cm, single packed	30 pieces

Sim-Plex

One-piece tracheostomy tube holder

Product properties

- 1,5 cm wide tracheostomy tube holder
- Strap length adjustable from 26-46 cm
- Hook fastener
- Easy and quick to use
- Universal - due to infinitely variable adjustment suitable for different neck circumferences
- Skin-friendly, soft materials
- Very thin and soft strap for discreet fixation
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Elastic strap
- Latex-free



REF	Article	Description	SU
200 446	Sim-Plex	Length range 26-46 cm, width 1,5 cm	60 pieces

Duo-Clip®

Two-part tracheostomy tube holder

Product properties

- Tracheostomy tube holder with soft 2,5 cm wide padding
- Strap length variably adjustable
- Available with hook or hook and loop fastener
- Easy and quick to use
- Suitable for all neck circumferences
- Skin-friendly, soft materials
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft foam with skin-friendly lamination
- Latex-free



REF	Article	Description	SU
100 801	Duo-Clip® H	Length range 37- 46 cm, 32,5 + 10 cm padding, hook fastener, single packed	30 pieces
200 801	Duo-Clip®H	Length range 37- 46 cm, 32,5 + 10 cm padding, hook fastener	60 pieces
100 802	Duo-Clip® K	Length range 33 - 42,5 cm, 32,5 + 10 cm padding, hook and loop fastener, single packed	30 pieces
200 802	Duo-Clip® K	Length range 33 - 42,5 cm, 32,5 + 10 cm padding, hook and loop fastener,	60 pieces

Duo-KidsClip

Two-part tracheostomy tube holder for children

Product properties

- Tracheostomy tube holder with soft 1,8 cm wide padding
- Strap length variably adjustable from 20-29,5 cm
- Available with hook and loop fastener
- Easy and quick to use
- Suitable for all neck circumferences
- Skin-friendly, soft materials
- Suitable for MRI
- Disposable product

Patient target group

- Neonates
- Infants
- Small children
- Children

Material

- Soft foam with skin-friendly lamination
- Latex-free



REF	Article	Description	SU
200 805	Duo-KidsClip	Length range 20 - 29,5 cm, 19,5 + 10 cm padding, hook and loop fastener	60 pieces
100 805	Duo-KidsClip	Length range 20 - 29,5 cm, 19,5 + 10 cm padding, hook and loop fastener, single packed	30 pieces

Duo-Fix

Two-part tracheostomy tube holder

Product properties

- Tracheostomy tube holder with soft 3 cm wide padding
- Strap length variably adjustable from 35-46 cm
- Easy and quick to use
- Skin-friendly, soft materials
- Hook and loop fastener with rounded hook and loop strap for quick and easy closing
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft foam with skin-friendly lamination
- Latex-free



REF	Article	Description	SU
200 803	Duo-Fix	Length range 35 - 46 cm, 25 + 10,5 cm padding	60 pieces

NOVO-Pad

Tracheostoma dressing

Product properties

- Tracheostoma dressing with zig-zag slit for easy change and secure hold with cannula in place
- Available in 2 sizes and 2 different thicknesses
- Easy to use
- Suitable for all common cannulas
- Especially recommended for just operated, irritated and exudating tracheostomata. Suitable for all stomas.
- No encrustation or adhesion to the wound surface
- Skin-friendly, soft materials, breathable
- High absorption capacity
- Suitable for MRI
- Single packed / Disposable product
- Sterile

Patient target group

- Adolescents
- Adults

Material

- Polyurethane film coating, polyurethane foam
- Latex-free



REF	Article	Description	SU
70-10	Novo-Pad	Width 9,5 cm, length 8,2 cm, thickness 0,6 cm, individually sterile packed	25 pieces
70-30	Novo-Pad 2	Width 14 cm, length 8 cm, thickness 0,4 cm, individually sterile packed	25 pieces

NOVO-Kids Pad

Tracheostoma dressing for children

Product properties


- Tracheostoma dressing with zig-zag slit for easy change and secure hold with cannula in place
- Easy to use
- Suitable for all common cannulas
- Especially recommended for just operated, irritated and exudating tracheostomata. Suitable for all stomas.
- No encrustation or adhesion to the wound surface
- Skin-friendly, soft materials, breathable
- High absorption capacity
- Suitable for MRI
- Single packed / Disposable product
- Sterile

Patient target group

- Small children
- Children
- Adolescents

Material

- Polyurethane film coating, polyurethane foam
- Latex-free



REF	Article	Description	SU
70-20	Novo-Kids Pad	Width 6,5 cm, length 6,5 cm, thickness 0,6 cm, individually sterile packed	25 pieces

NOVO SealPad ○

Self-adhesive and sterile wound dressing

Product properties

- Wound localization device (blue button) for optimal pressure when talking, coughing, sneezing as well as during phonation and breathing exercises
- Optimal adaptation to the anatomical conditions of the body
- Single packed / Disposable product
- Sterile

Patient target group

- Children
- Adolescents
- Adults

Material

- Silicone wound contact layer, polyurethane foam pad
- Waterproof PU film backing



REF	Article	Description	SU
TS003	Novo Seal Pad, Oval	Dimensions 7,5 x 6,5 cm	10 pieces
TS005	Novo Seal Pad, Square	Dimensions 7,5 x 7,5 cm	10 pieces



1.5 Inhalation & Aerosol Therapy Nebulizers



Nebulizer Kit ○

With mouthpiece

Product properties

- Kit consists of nebulizer with mouthpiece, T-piece, 15 cm corrugated tube and 2,1 m O₂ tube
- Single packed / Disposable product

Patient target group

- Children
- Adults

Material

- PVC



Nebulizer Pot	6 ml volume
O ₂ Tube	210 cm with push-in and screw cap
Corrugated Tube	as reservoir
Flow Rate (gpm)	0.40±0.05@7LPM
Min. Flow Rate	1.85±0.30LPM
Droplet Size (MMAD)	3,30 μ ±1,00@7LPM

REF	Article	SU
IC-VMU 01	Nebulizer kit with mouthpiece	100 pieces

Nebulizer Kit ○

With mouthpiece

Product properties

- Kit consists of nebulizer with mouthpiece, T-piece with non-return valve and norm connectors and 2 m O₂ tubing
- Single packed / Disposable product

Patient target group

- Children
- Adults

Material

- PVC



Nebulizer Pot	6 ml Volume
O ₂ Tube	200 cm with push-in and screw cap
T-piece	also suitable for tracheal cannula

REF	Article	SU
30148-WOC	Nebulizer kit with mouthpiece	50 pieces
30148-V12.	Nebulizer kit, T-piece with non-return valve, cover cap	50 pieces
30512	Nebulizer pot	200 pieces

Nebulizer Kit ○

With mask and O₂ tube

Product properties

- Kit consists of nebulizer with nebulizer pot, mask and 2,10 m O₂ tube
- Single packed / Disposable product

Patient target group

- Neonates
- Small children
- Children
- Adolescents
- Adults


Material

- PVC



Nebulizer Pot	6 ml Volume
O ₂ Tube	210 cm with plug-in and screw cap 9/16'-9/16' thread & 6 mm funnel neck
Nebulizer	18 mm OD & 6 mm funnel attachment
Mask	18 mm ID
Flow Rate (gpm)	0.40±0.05@7LPM
Min. Flow Rate	1.85±0.30LPM
Droplet Size (MMAD)	3,30 μ±1,00@7LPM

REF	Article	SU
IC-VMA 03	Nebulizer kit for babies	100 pieces
IC-VMA 01	Nebulizer kit for children	100 pieces
IC-VMA 02	Nebulizer kit for adults	100 pieces



1.6 Patient Positioning

Head positioning device



NOVO PO QView

Head positioning device with mirror

Product properties

- Head positioning device with mirror for patient observation
- Shatterproof foil mirror
- Adhesive strips enclosed for attachment to the mirror
- Optionally available in 12 cm or 15 cm height
- Skin-friendly medical tape for forehead and chin support
- Dimensions (L x W x H) 30 x 24 x 12 or 30 x 24 x 15 cm
- Disposable product/single packed

Patient target group

- Adults

Material

- Polyurethane (PU), polyacrylic
- Latex-free

REF	Article	SU
Q3VIEW	NOVO PO QView Head Positioning Devices 12 cm	8 pieces
QPHPD-ADULT1	NOVO PO QView Head Positioning Devices 15 cm	6 pieces



NEW



2. Patient Safety

On the following pages you will find products related to the topic of fixation and pressure relief, including fixations for patient self-protection as well as cushionings for tubing systems.

2.1 Decubitus Protection	75
2.2 Phototherapy	76
2.3 Patient Self-Protection	77

NOVO PFT 100

Cushioning for tubing systems

Product properties

- For tubing systems of various sizes
- Padding length 7 cm, outer diameter 1 cm
- Easy and quick to use
- Can be shortened individually
- Protection against decubitus ulcer
- Suitable for MRI
- Disposable product

Patient target group

- For patients of all ages

Material

- Foam
- Latex-free



REF	Article	Description	SU
40-10	NOVO PFT 100	Length 7 cm, outer diameter 1 cm, adhesive strip enclosed	120 pieces

Ute-Fix

Cushioning for probes and drainages

Product properties

- Cushioning with adhesive strip
- For probes, drainages and tubing systems of various sizes
- Available in 3 sizes
- Easy and quick to use
- Universally applicable due to different sizes
- Light and soft material
- Suitable for MRI
- Disposable product

Patient target group

- For patients of all ages

Material

- Foam
- Latex-free



REF	Article	Description	SU
40-06	Ute-Fix 1	Length 2 cm / OD 1,5 cm	60 pieces
40-05	Ute-Fix 2	Length 2 cm / OD 2 cm	60 pieces
40-01	Ute-Fix 3	Length 2 cm / OD 2,5 cm	60 pieces

Eye-Care Kids



One-piece eye protection

Product properties

- One-piece eye protection
- Available in 2 sizes
- Hook and loop fastener
- Easy and quick to use
- Skin-friendly and soft materials
- Light weight
- Suitable for MRI
- Disposable product



Patient target group

- Neonates


Material

- Soft foam with skin-friendly lamination
- Latex-free

REF	Article	Description	SU
60-31	Eye-Care Kids 32	Length 32 cm	15 pieces
60-32	Eye-Care Kids 38	Length 38 cm	15 pieces

2.3 Patient Self-Protection

Hand and Foot Fixations



Hand-Fix 22/25/30



NEW

Hand and foot fixation

Product properties

- Soft padded cuff with fixation strap
- Available in various sizes and different versions
- Easy and quick to use
- Skin-friendly materials
- Secure hold thanks to strong hook and loop fastener and tear-resistant fixation strap
- Soft and thick padding to protect against decubiti and injuries
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

Material

- Soft foam with skin-friendly lamination, polyester
- Plastic (fixation ring)
- Latex-free



REF	Article	Description	SU
60-61	Hand-Fix 22	With fixation ring and fixation strap, 7 x 22 cm (w x l), strap length 150 cm (enclosed)	16 pieces
60-63	Hand-Fix 25	With fixation ring and fixation strap, 7 x 25 cm (w x l), strap length 150 cm, (enclosed)	16 pieces
60-65	Hand-Fix 30	With fixation ring and fixation strap, 7 x 30 cm (w x l), strap length 150 cm, (enclosed)	16 pieces

Hand-Fix Kids 12/17



Hand and foot fixation for children

Product properties

- Soft padded cuff with fixation strap
- Available in various sizes and different versions
- Easy and quick to use
- Skin-friendly materials
- Secure hold thanks to strong hook and loop fastener and tear-resistant fixation strap
- Soft and thick padding to protect against decubiti and injuries
- Suitable for MRI
- Disposable product

Patient target group

- Children

Material

- Soft foam with skin-friendly lamination, polyester, plastic
- Latex-free



REF	Article	Description	SU
60-01	Hand-Fix Kids 12	With sewn-on fixation strap, 3 x 12 cm (w x l), strap length 2x65 cm, for neonates, infants and small children	30 pieces
60-08	Hand-Fix Kids 17	With sewn-on fixing strap and eyelet, 4 x 17 cm (w x l), strap length 2x69 cm, for small children and children	30 pieces
60-50	Hand-Fix Kids Mix	14 x Hand-Fix Kids 12, 14 x Hand-Fix Kids 17, for neonates, infants, small children and children	28 pieces



3. Infusion and Surgical Drain Management

Optimally fitting fixations support the medical devices in fulfilling their tasks and are thus the basis for safe infusion & drainage management. Here you will find fixation solutions for infusions, drains and decubitus ulcer prophylaxis.

3.1 Fixations

82

ECMO-Fix

Fixation for catheter tubing systems

Product properties ECMO-Fix 1

- Patented product solution
- Fixation consisting of 1 fixation strap, 2 fixation plasters and 1 fixing loop
- Fixation strap with 2 hook and loop straps and hook and loop fastener
- Easy to use
- Fixation strap individually adjustable to head circumference
- Position correction of the catheter tubing set possible at any time
- Fixing loop alignable to the position of the acces
- Optimal cushioning of the catheter tube from the skin
- Disposable product
- Suitable for MRI



Product properties ECMO-Fix 2

- Fixation consisting of 1 fixation strap and 1 plaster fixation
- Fixation strap with 2 hook and loop straps and hook and loop fastener
- Easy to use
- Fixation strap individually adjustable to the thigh circumference
- Position correction of the catheter tubing set possible at any time
- Skin-friendly, soft materials
- Optimal cushioning of the catheter tube from the skin
- Extra wide strap to protect against decubitus ulcer
- Disposable product
- Suitable for MRI




Patient target group

- Adolescents
- Adults


Material

- Soft foam with skin-friendly lamination
- Latex-free



ECMO-Fix 1


For use during veno-venous ECMO via the jugular vein in the head area



ECMO-Fix 2

For use via the femoral vein on the leg and for iLA/PECLA procedures

REF	Article	Description	SU
100-300-01	ECMO-Fix 1	Width 2,5 cm, length 72 cm, thickness 0,5 cm	10 pieces
100-300-02	ECMO-Fix 2	Width 4 cm, length 72 cm, thickness 0,5 cm	10 pieces



Main-Lock® 1-3 / Main-Lock® 9

Plaster fixation

Product properties


- Available in different sizes
- Easy and quick to use
- Particularly stable hook and loop fastener
- Skin-friendly materials
- For strain relief at the puncture site
- Plaster on the surface hydrophilic
- Suitable for MRI
- Single packed / Disposable product

Patient target group

- For patients of all ages

Material

- Foam, acrylic adhesive layer
- Non-woven material, acrylic adhesive layer (Main-Lock® 9)
- Latex-free



REF	Article	Description	SU
90-60-81	Main-Lock® 1	Length 7 cm, width 3,7 cm, hook and loop strap with interior adhesive surface	100 pieces
90-60-82	Main-Lock® 2	Length 9 cm, width 4,3 cm, hook and loop strap with interior adhesive surface	100 pieces
90-60-83	Main-Lock® 3	Length 9 cm, width 4,3 cm, plaster with 2 separate hook and loop straps with interior adhesive surface	100 pieces
90-50-09	Main-Lock® 9	Length 9 cm, width 4 cm, hook and loop strap with interior adhesive surface	100 pieces

Main-Lock® Set

Plaster fixation for catheters and drainages

Product properties

- Fixation consisting of a plaster and a tube fixation
- Can be removed from tubes without leaving any residue
- Multiple fixations and change of tubing course possible
- Also suitable for use with long-term catheters and drainages
- For strain relief at the puncture site
- Easy and quick to use
- Skin-friendly materials
- Suitable for MRI
- Single packed / Disposable product

Patient target group

- Small children
- Children
- Adolescents
- Adults

Material

- Non-woven material, acrylic adhesive layer
- Latex-free



REF	Article	Description	SU
90-30-15	Main-Lock® Set 10	Length 10 cm, width 4 cm, for small children, children, adolescents and adults	40 pieces
90-30-20	Main-Lock® Set 14	Length 13,5 cm, width 4 cm, for children, adolescents and adults	40 pieces

Nasal-Fix

Nasal fixation for probes

Product properties

- Suitable for all common tube sizes
- Easy and quick to use
- Thin and flexibly adaptable material
- Suitable for MRI
- Single packed / Disposable product

Patient target group

- Adolescents
- Adults

Material

- Non-woven material, acrylic adhesive layer
- Latex-free



REF	Article	Description	SU
100-200-01	Nasal-Fix	Length 8,2 cm, width 7 cm	80 pieces

Splint-Fix

Stabilisation and positioning splint

Product properties

- Skin-friendly padded securing splint with flexible aluminium core
- Sewn on fixation straps with hook and loop fastener
- Easy and quick to use
- Easily adjustable splint
- Different sizes allow an individual use
- MR conditional
- Disposable product

Patient target group

- For patients of all ages

Material

- Soft foam with a skin-friendly cover and aluminium core, soft rubber
- Latex-free



Various application possibilities



Splint-Fix 9

To ensure a secure stabilisation the straps are used crosswise.



Splint-Fix 13

Immobilises the wrist and allows a clear view of the puncture site at the same time.



Splint-Fix 22

Protects the access when mobilising and repositioning the patient.

REF	Article	Description	SU
60-21	Splint-Fix 9	Tapered splint: 9 x 3-4 cm (l x w), 2 fixation straps: 14 x 1.6 cm (l x w), for neonates, infants and small children	20 pieces
60-22	Splint-Fix 13	Splint: 13 x 4 cm (l x w), 2 fixation straps: 14, 16 x 2.5 cm (l x w), for neonates, infants and small children	20 pieces
60-23	Splint-Fix 16	Splint: 16 x 5 cm (l x w), 2 fixation straps: 17, 19 x 2.5 cm (l x w), for infants and children	20 pieces
60-11	Splint-Fix 22	Splint: 22 x 9 cm (l x w), 2 fixation straps: 30, 40, x 5 cm (l x w), for children, adolescents and adults	14 pieces
60-13	Splint-Fix 33	Splint: 33 x 8,5 cm (l x w), 3 fixation straps length 25, 32 and 42 x 5 cm (l x w), for children, adolescents and adults	10 pieces



Clip-Fix

Fixation for cables and tubes

Product properties

- Easy to use
- Universally applicable
- Strong hook and loop
- Clip-Fix 2 Trio with 3 Velcro straps for fixing several tubes and cables
- Clip-Fix 1 Duo with clip
- 2.5 cm wide hook and loop strap, with one or two steel clips depending on the version


Clip-Fix is reusable. Reprocessing by washing at 60 °C with a commercially available detergent, washable up to 10 x.

Patient target group

- For patients of all ages

Material

- Plastic, stainless steel
- Latex-free



REF	Article	Description	SU
30-30-01	Clip-Fix 1	Width 2,5 cm, length 11 cm, 1 steel clip	50 pieces
30-30-02	Clip-Fix 2	width 2,5 cm, length 17,5 cm, 1 steel clip	50 pieces
30-30-03	Clip-Fix 3	width 2,5 cm, length 24 cm, 1 steel clip	40 pieces
30-30-04	Clip-Fix 1 Duo	width 2,5 cm, length 13 cm, 2 steel clips	50 pieces
30-30-08	Clip-Fix 2 Trio	Width 2,5 cm, length 17 cm, 3 hook and loop straps	50 pieces



Clip-Fix 1
Guiding and holding of several tube systems with fixation at the bed sheet.



Clip-Fix 3
Guiding and holding of the breathing circuit and the cable systems at the service rail placed on the handgrip of the rescue vehicle.

Uni-Fix

Fixation for cables and tubes

Product properties

- One-piece fleece strap with hook and loop fastener and eyelet
- Available in 2 sizes
- Easy and flexible to use
- Adjustable in length
- To reduce tensile stress
- Suitable for MRI

Uni-Fix is reusable. Reprocessing by washing at 60 °C with a commercially available detergent, washable up to 10 x.

Patient target group

- For patients of all ages

Material

- Plastic, stainless steel
- Latex-free



REF	Article	Description	SU
30-50-01	Uni-Fix 1	Length 24,5 cm, width 2 cm	20 pieces
30-50-02	Uni-Fix 2	Length 40,5 cm, width 2 cm	20 pieces



4. Monitoring

For a safe patient care, we also offer products in the field of patient monitoring. On the following pages you will find, for example, fixations for transducers or sensors as well as measuring devices.

4.1 Blood Pressure Measurement	92
4.2 Cuff Pressure Measurement	93
4.3 SPO ₂ Measurement	94

Transducer-Fix



Fixation for transducers

Product properties

- Cuff with 2 hook and loop straps to fix a maximum of 2 measuring cables
- Available in 2 versions
- Easy to use
- Ideal for transport or mobile patients
- Extra wide for stable fixation
- Soft padding
- More mobility because transducers remain securely fixed when moving
- Suitable for MRI
- Disposable product

Patient target group


- Adolescents
- Adults

Material

- Soft foam with skin-friendly lamination
- Latex-free



REF	Article	Description	SU
30-45	Transducer-Fix H	Length 56 cm, width 12 cm, 2 horizontal hook and loop straps	36 pieces
30-45-03	Transducer-Fix V	Length 56 cm, width 12 cm, 2 vertical hook and loop straps	36 pieces



Pylant Monitor



Cuff pressure manometer


Product properties

- Monitor for measuring and checking the cuff pressure of endotracheal tubes and tracheal cannulae
- Luer-Slip male
- Suitable for MRI
- Disposable product

Patient target group

- Adolescents
- Adults

REF	Article	Description	SU
20-60	Pylant Monitor	Cuff pressure manometer	12 pieces



PulsOxy-Fix Kids



Fixation for pulse oximeter sensors

Product properties

- One-piece pulse oximeter sensor fixation
- Available in different lengths and widths
- Easy and quick to use
- Fixation on fingers, feet and other limbs possible
- Additional internal adhesive strip to secure the sensor with PulsOxy-Fix Kids 12 and 17
- Skin-friendly materials
- Non-adhesive attachment to the skin - gentle on the sensitive skin of the little patients
- Suitable for MRI
- Disposable product





Patient target group

- Premature babies
- Neonates
- Infants
- Small children
- Children

Material

- Soft foam with skin-friendly lamination
- Latex-free

REF	Article	Description	SU
30-14	PulsOxy-Fix Kids 10	Length 10 cm, width 1.8 cm, for premature babies, neonates and infants	120 pieces
30-10	PulsOxy-Fix Kids 12	Length 12 cm, width 3 cm, interior adhesive strip for sensor protection, for neonates, infants and small children	80 pieces 
30-12	PulsOxy-Fix Kids 17	Length 17 cm, width 4 cm, interior adhesive strip for sensor protection, for small children and children	60 pieces 





5. Hygiene & Cleansing

Oral hygiene is an important factor in the treatment of intensive care patients. Here we present our innovative product solution for the care and protection of the oral flora.

5.1 Oral/Dental care

98

Hygiene & Cleansing

Gentle and effective product solutions for daily care and cleaning

Our Jack Pro® Products with the Oral Hygiene Concept® enable carefully performed antiseptic oral care. In accordance with the Commission Recommendation (RKI), we use substances with proven efficacy.

Jack Pro® Liquid plus and **Gel plus** are bioscientific products for oral cleansing, germ reduction and wound care in respiratory patients, as well as in oncological, paediatric and geriatric care.

Our product solutions support both the self-management competence of patients and caring relatives and the care process for the maintenance, promotion and treatment of oral health by professional nursing staff in inpatient and outpatient facilities. This includes the following target groups in particular:

- Laryngectomy or tracheostomy patients
- Elderly patients, adults and children in need of care, including patients in intensive care, oncology, radiology and oral and maxillofacial surgery
- Patients with immune deficiency
- Patients with periodontitis-associated diseases

The main pathogens of classic nosocomial infections are combated. They are already inactivated, weakened in growth or killed with a low dose of active ingredient. Jack Pro® protects and supports the physiological oral flora.

Product properties

- Mild ingredients - gentle on the oral cavity, teeth and gums
- Long-lasting use - suitable for daily use
- Biocompatible - maintains the ecological balance of the oral cavity and supports its regeneration
- The gel can be used to clean the teeth, it does not need to be be rinsed out
- Creates a long-lasting moist oral mucosa
- No protein defect in bleeding and purulent wound healing disorders
- As an oral cleansing concept, it complies with the bundle recommendation of the Cambridge VAP Guideline of May 2022

Composition

1. Rinse:

0,025 g Polyaminopropyl Biguanide; 0,005 g Tosylchloramide Sodium; Phenylalanine; Xylitol; Ethylenediaminetetraacetate; Allantoin; Aroma (Xylitol); C.I. 42051"

2. Gel:

0,33g Sodium Fluoride; 0,1g Extract. Salic. Alba Fluid; Xanthan; Aloe Vera; Glutamic Acid; NN-DiaceticAcid Tetrasodium Salt; Aroma; Menthol



REF	Article	Description	SU
40020	Jack Pro® Gel plus	Bottle 100 ml	12 pieces
40030	Jack Pro® Liquid plus	Bottle 250 ml	6 pieces
41005	Jack Pro® 24-36-48h UNIT	2x Jack Pro® Gel plus 30 ml 1x Jack Pro® Liquid plus 35 ml in dispenser box	30 pieces
42007	Swabby	Oral swab	100 pieces



Novo Klinik-Service GmbH
Technologiepark West – Haus 15
Zum Frenser Feld 1
50127 Bergheim • Germany

Tel. +49 2271 99404-0
Fax +49 2271 99404-700
info@novo-med.de
www.novo-med.de

